

Washing Buffer

**Washing Buffer for ELISA, EIA, RIA,
Protein Assays, Western Blotting,
Immuno-PCR and Immunohistochemistry**

Washing Buffer saves time and money!

Washing Buffer is useable for all kinds of Immunoassays.

Fields of application are ELISA, EIA, RIA, Protein Arrays, Western Blotting and Immuno-PCR.

For Immunohistochemistry also available without Tween.

Useable with ELISA washers



Washing Buffer

Washing steps are needed to remove unbound and excessive components, which are able to interfere with the assay.

Useable with automatic washers.

(Please follow instruction manual of the washer!)

For Immunohistochemistry we offer a Washing Buffer without Tween and detergents.

10x stock solution

dilute with water!

composition of
work dilution:

Buffer based on PBS with Tween and 350 mM NaCl

pH-value:

pH 7,4

preservative:

contains 0,005% thimerosal

storage:

shelf life at 2-8 °C, 6 months

use working dilution immediately

For research use only, not for diagnostic use

available

packaging size:

Washing Buffer

500 ml order number 140500

Washing Buffer without Tween for Immunohistochemistry

500ml order number 141500